

Disposable Microforceps: Eckardt End-gripping

(25 gauge / 0.5 mm)

Article No: 1286.WD05



Mr Danny Mitry
Consultant Ophthalmologist,
UK, Royal Free London
NHS Foundation Trust

“The 25G DORC Eckardt Forceps are the most versatile forceps I have come across. I find them ideal in all types of macula surgery from ILM peel to adherent PVR membranes.”

- **Microforceps with ‘nail-like’ end-gripping platforms**
- **For precise grasping and improved visualization of thin membranes in the macular region**
- **Developed with Prof. M.D. Claus Eckardt**

[Visit product page](#)

