


Disposable Eckardt TwinLight Chandelier (27 gauge / 0.4 mm)

Ref. 3269.MBD27



Umberto Lorenzi MD, FEBO
Head of Medical and Surgical Retina
Hospital Center University
Charles-Nicolles Rouen, France

“To achieve the best surgical result we need excellent visualization.

In vitreoretinal surgery hands-free light allows for better control at all times, accurate shaving of the vitreous base whilst indenting through 360°, and using bimanual approach to solve complex cases.

Disposable Eckardt TwinLight Chandelier 27G is a very good light fiber that doesn't need specific trocars and can be inserted through the conjunctiva using a guidance needle.

EVA LED light source in conjunction with Eckardt TwinLight Chandelier provides comfortable visualization to perform all vitreoretinal cases.”

