


DISPOSABLE I/A HANDPIECES WITH TEXTURED GRIP

Coaxial and bimanual



DISPOSABLE I/A HANDPIECES WITH TEXTURED GRIP

- Textured handpiece surface for comfortable handling / grip
- Compatible with standard luer connector
- Colour coding for each Gauge size



Disposable Coaxial I/A handpieces



23G Coaxial

Compatible with EquiPhaco 1.8mm



21G Coaxial

Compatible with EquiPhaco 2.2mm



21G Coaxial Soft Tip

Compatible with EquiPhaco 2.2mm



20G Coaxial

Compatible with EquiPhaco 2.4mm



- Soft contact point to fragile tissue and IOL
- Soft tip gives optic feedback to the surgeon
- Soft aspiration port

Disposable Bimanual I/A handpieces



0.8mm



1.1mm



Coaxial

PRODUCT CODE	SIZE	INCISION SIZE	ASPIRATION CANNULA	ASPIRATION PORT	ASPIRATION HEAD
SMS370P	23G / 0.65mm	1.8mm	angled 45°	0.30mm	polished
SMS355P	21G / 0.8mm	2.2mm	angled 45°	0.35mm	polished
SMS390S	21G / 0.8mm	2.2mm	angled 45°	0.40mm	silicone tip (soft-tip)
SMS325P	20G / 0.9mm	2.5mm	angled 45°	0.35mm	polished

Bimanual

PRODUCT CODE	SIZE	INCISION SIZE	ASP. CANNULA	ASP. PORT	ASP. HEAD	IRR. CANNULA	IRR. PORT
SBS320A	23G / 0.65mm	0.8mm	aeroplane shape	0.30mm	textured	aeroplane shape	2 ports, 0.35mm
SBS310A	21G / 0.8mm	1.1mm	aeroplane shape	0.35mm	textured	aeroplane shape	2 ports, 0.45mm