



**Disposable Microforceps:**  
**Extended reach**  
**Wide grip, 27 gauge / 0.4 mm**

Ref. 1286.WRD04



**Developed in collaboration with**  
**Prof. Fanis Pavlidis, Germany**

“The surgical comfort and efficiency of the 27G peelings are mostly dependant on the ergonomics and stiffness of the forceps. The Extended Reach - Wide Grip 27G forceps combines significantly increased rigidity for direct response and a wider gripping tip allowing for a reduced number of grasp actions and minimized risk of retinal trauma. In addition, the stronger end gripping tip with serrations allows for stronger grip of membranes.”

**Features & Benefits in comparison**  
**with the older 27G forceps:**

- 60% increased gripping width
- 260% improvement in gripping area
- 40% enhanced stiffness