



SIL-2000-S

Silicone Oil

Ref. SIL-2000-S



Prof. Anne-Catherine Gribomont,
Vitreoretinal surgeon at UCL St Luc,
Brussels, Belgium

“Since DORC launched their Silicone oil 2000, I have switched from the 5000 cSt silicone oil to the 2000 cSt. The main reason for switching is the fact that both - injection and extraction - are much easier and faster with the Silicone oil 2000. Also, when using the cannula system and silicone oil you always have some oil that sticks to the metal of the cannula. Although this is also the case with the Silicone oil 2000, the removal is much easier. For me the efficiency of injecting and removing of the 2000 cSt oil is obvious.”

Features & Benefits:

- Single compound oil
- Narrow viscosity range
- Highly purified
- Low levels of emulsification
- Available in 10ml syringe
- Faster to infuse/extrude than 5000cSt oil