



# Disposable Vitrectomy Lenses

Article No: 1284.DD/GD/ID



**Professor Dr.**  
**Angeles Hernandez Vazquez,**  
Nuestra Señora de la Luz  
Mexico City

**“Vitrectomy lenses provide me with enough magnification to clearly see segmentation planes or peel membranes with a clear view on details, with perfect stereoscopy for precise movements; it is a very useful tool for delicate or complicated cases.”**

- Disposable Vitrectomy Lens: Flat (1284.DD):  
- view of the central vitreous and fundus



[Watch surgical video](#)

- Disposable Vitrectomy Lens: Wide View (1284.GD):  
- view from the central up to the equatorial portion of the vitreous cavity



[Watch surgical video](#)

- Disposable Vitrectomy Lens: 30° Prism (1284.ID):  
- view of mid (20°) or far (30°) periphery



[Watch surgical video](#)